


C A G E N I X™

ACCUFRA[®] IC™

CUSTOM CERAMIC CROWNS ON OUR ACCUFRA[®] FRAMEWORK FOR ABOUT THE SAME PRICE AS A TRADITIONAL FIXED-HYBRID.



HOW DOES ACCUFRA[®] IC STACK UP

ADVANTAGES OVER OTHER COMMON RESTORATIVE OPTIONS

VS. FIXED-HYBRIDS

SUPERIOR AESTHETICS

Full-contour crowns available in zirconia with optional cutbacks for porcelain stacking

SUPERIOR DURABILITY

Custom crowns perfectly mated to our award-winning AccuFrame Plus Titanium frameworks

SUPERIOR FUNCTION

Patented crown-on-post design transfers occlusal loads, improving chewing performance and gingival durability

SUPERIOR REPAIRABILITY

Quick and easy crown replacement via our lingual retrieval notch and digitally stored crown files

VS. INDIVIDUAL CROWNS

SUPERIOR VALUE

AccuFrame IC pricing offers substantial savings over traditional individual crown restorations

SUPERIOR PROCESSING

Easy system for processing acrylic gingiva means savings both in time and money

SUPERIOR TURN-AROUND

Fewer steps and patient visits plus quicker manufacturing turnaround means lower costs and quicker delivery

SUPERIOR REPAIRABILITY

Our patented crown-on-post design significantly improves ability to remove/repair the restoration

ACCUFRA[®]ME IC

INDIVIDUAL CROWN SYSTEM

Introducing AccuFrame IC

For years, the only options for fixed-hybrid cases have been traditional restorations with common denture teeth that are affordable but lack durability or extremely expensive, complex, individual crown options. Our AccuFrame IC (Individual Crown) System is removing the price barrier to improved look and performance of your fixed-hybrid cases — and making it possible for you to offer the highest quality full-arch restorations at truly affordable prices.

AccuFrame IC restorations utilize full-contour individual crowns in high-strength zirconia or layered porcelain, on our award-winning AccuFrame Plus titanium frameworks. The result? Guaranteed fit, superior appearance and long-lasting performance.



AccuFrame IC vs. Other Options

As you can see below, other restorations lack a means for the denture tooth forces to directly engage the bone structure. AccuFrame IC's titanium pontic form and framework structure has been specifically engineered to mate perfectly with the custom crowns, thereby transferring the occlusal force directly to the bone. This results in the best-performing, longest-lasting restoration available.

What type of patients benefit from Individual Crowns

- Patients seeking most realistic appearance available
- Patients seeking to minimize return visits for repair/rework
- Patients seeking to regain full functionality of all teeth
- Patients with increased wear traits (like Bruxers or grinders)
- Patients wanting the highest quality restoration for their money



Traditional Denture
No structural contact between occlusal surface and bone



Fixed-Hybrid
Frame engages bone but teeth supported via acrylic gingiva only



AccuFrame IC
Direct transfer of occlusal loads to bone for optimal performance



For more info call (866) 964-5736 or visit cagenix.com