

ACCUFAME 360

A Bar For Life

PMMA

- AccuFrame 360 Titanium Framework
- PMMA monolithic overlay structure
- Acrylic gingival contours

Perfect alternative to fixed-hybrids for price-sensitive or high-wear patients. Spend less time and money on repairing traditional fixed-hybrids.

Universal Framework

Custom milled out of 6AL-4V Titanium with all current or future overlay options engineered to be easily interchangeable.

LIFETIME
FRAMEWORK WARRANTY

Future Proof

Why lock your patients into today's restorations when tomorrow's advances are just around the corner?

- Our material-neutral system offers the widest range of current and future restorative options.
- Universal frame interface allows for easy future upgradability.

Zirconia

- AccuFrame 360 Titanium framework
- Zirconia monolithic overlay structure
- Layered pink porcelain gingiva
- Optional porcelain on buccal edge

Premium durability and aesthetics for patients who demand the best. An easy upgrade from other system options.




One Frame. Endless Possibilities.

One Frame. Endless Possibilities.

ACCUFRAME 360

A Bar For Life



Monolithic PMMA Overlay


Frame: 6AL-4V Titanium

Tooth: PMMA Surface

Overlay: PMMA Monolithic

Gingiva: Acrylic

An ideal alternative to fixed-hybrids for those seeking durable or easily replaceable restorations. This monolithic PMMA lends itself for use in bruxers or for value-conscious patients. The interchangeable overlays streamline repairs or upgrades.



Monolithic Zirconia Overlay

Frame: 6AL-4V Titanium

Tooth: Zirconia or Porcelain

Overlay: Zirconia Monolithic

Gingiva: Porcelain

Upgrade from your old monolithic restorations with our improved durability, simple retrievability and 30% lower cost. Porcelain stacking on gingival and optional buccal surfaces offers best-in-class aesthetics at an incredible value.

Additional Benefits:

- As little as 8mm restorative space required
- Lifetime framework warranty
- Choose from full wrap or metal-base options
- Support for any millable material available
- Innovative patented design

C A G E N I X™

Visit cagenix.com or call (866)964-5736 for more information.

v0621.01